

DCM Documentation

Requirements

- The latest version of Python
- Pillow
 - Install with `pip install pillow`

Instructions for launching DCM

1. Clone the repository
2. Open terminal and navigate to the DCM_group44 folder.
`cd DCM_group44`
3. Launch the DCM by running main.py through the terminal.
`python3 main.py`

Explanation of each file and class organization system

Main.py (window): the start file or “welcome page” that can call the login window or registration window upon the respective button being clicked.

Login.py (window): launched from the main file. Checks if inputted credentials are valid and exist in the database, and launches modeSelection, passing on the username so that the modeSelection window knows what row in the database to import values from and save values to.

Registration.py (window): launched from the main file. Checks if inputted credentials are valid and saves them to the database. It launches modeSelection, passing on the username for the same reasons as Login.py above.

Patient.py: defines the Patient class whose instance stores the values and parameters for a patient, as well as methods to copy values from the database and save values to the database.

data.py: defines the createDB() function which is called by several window scrips to initialize a database file if it doesn't already exist.



modeSelection.py (window): launched from either registration.py or login.py. Instantiates a Patient object with default parameter values and calls the method Patient.copyFromDB() which overwrites those default values with any custom values in the database that relate to the patient with that username. The window shows 4 modes to choose from: AOO, VOO, AAI, VVI. When one of the modes is chosen and “Next” button is clicked, it launches a specific frame from the 4 frames in pacingModes.py. Transitions between modeSelection.py and pacingModes.py are switching frames

instead of switching windows. When launching `pacingModes`, the `Patient` object stored in this class is transferred to the new pacing mode frame launched, giving `pacingModes.py` the info needed about the patient in order to copy from and to save to the database.

`pacingModes.py (frames)` : launched from `modeSelection.py` from which it gets a `Patient` object with the values relating to a specific patient in the database. Four mode frames exist: `AOO`, `VOO`, `AAI`, `VVI`, which are all children of a parent class `PacingMode`. The parent class `PacingMode` defines methods to add each label and entry, and the children objects determine which of those methods to call depending on which parameter entries are relevant to that mode. Each frame displays the specific fields, labels, and checkboxes required for the mode it represents. Upon launching, it copies all the values from the `Patient` object received from the `modeSelection` window to the entries in that frame. In other words, when a user logs in, they would find their own values in the entries instead of the default values. When the “Confirm” button is clicked, all entry values are transferred to a `Patient` object, and `Patient.saveToDB()` method is called, transferring the patient information to the database.

Demo

Some of our favourite pages!


<p>Welcome Page</p>	 <p>The screenshot shows a web browser window titled "Pacemaker". The status bar at the top indicates "Not connected - No new device". The main content area features the word "Welcome" in a large, bold, black font. Below the text is a large, stylized heart graphic. The heart is composed of two parts: a solid, textured red heart on the left and a red circuit board pattern on the right. At the bottom of the page, there are two buttons: "Login" and "Register".</p>
<p>Login Screen</p>	 <p>The screenshot shows a web browser window titled "Pacemaker Login". The status bar at the top indicates "Not connected - No new device". The main content area features the word "Login" in a large, bold, black font. Below the text is a smaller version of the stylized heart graphic. Underneath the heart, there are two input fields: "Username" and "Password". At the bottom of the page, there are two buttons: "Back" and "Login".</p>

Registration Page

Pacemaker Register

Not connected - No new device

Register



Username

Password

Re-enter password

Pacing Mode Selection Page

Pacemaker | Mode Selection

Not connected - No new device

Please choose a pacing mode

☒ AOO

☐ VOO

☐ AAI

☐ VVI

VVI Pacing Mode
Page: initially with
nominal values
specified in
“PACEMAKER”
document

Pacemaker | VVI Pacing Mode

Not connected - No new device

VVI

Please enter the values for the following parameters.
These values will be checked to ensure that they are valid entries.
If you have not set any values yet, they will be set to the nominal values.

Lower Rate Limit (ppm):	69
Upper Rate Limit (ppm):	120
Ventricular Pulse Width (ms):	0.4
Ventricular Amplitude (V):	3.5
Ventricular Sensitivity (mV):	2.5
Ventricular Refractory Period (ms):	320
<input type="checkbox"/> Hysteresis (ppm):	
Rate Smoothing (%):	0

[Back](#) [Confirm changes](#)

This is just an example. All of the other pacing modes look similar, just with their respective parameters and nominal values initially inputted.